


# ICE WRAP Cool & Hot Dual System



ICE cooling firmly wraps the RF heating energy generated in the deep layer and keeps the skin moderately comfortable during the RF high temperature. This cooling operation aids in maximizing the level of RF on the skin and decreases treatment time.

ICE wrap can help fat cells around the organs as well as under the dermis via the powerful deep layer RF energy. ICE wrap has three handles for use on any part of the body including lifting and wrinkle care for the face.




ICE cooling firmly wraps the RF heating energy generated in the deep layer and keeps the skin moderately comfortable during the RF high temperature.

| RET RF                    | 4-polar RF              | ICE cooling            |
|---------------------------|-------------------------|------------------------|
| - Frequency : 0.5 Mhz     | - Frequency: 0.5 Mhz    | - Frequency: Max.-20°C |
| - Handle size : 85π, 60π  | - Intensity: 1~10 level | - Intensity: 1~5 level |
| - Intensity: 1~10 level   | - Time: 1~99 mins       |                        |
| - Time: 1~99 mins         | (Basic: 20 mins)        |                        |
| (Basic: 20 mins)          |                         |                        |
| - Output power: 600W      |                         |                        |
| - Dimension: 300×405×1010 |                         |                        |
| - Weight: 35kg            |                         |                        |
| - Power consumption: 700W |                         |                        |



- Normal-big size by 80r
- Face tightening
- Skinny & thigh by 60r

# ICE WRAP Cool & Hot Dual System



A user can select a factory recommended mode according to a client's treatment part from MODE SELECTION, BASIC, LONG1, LONG2, SMALL1, SMALL2 and EASY(F) mode. Each mode has programmed by recommended RF and ICE level based on users' effective result and satisfactions. If a user selects a MANUAL mode, it can be programed manually with UP/DOWN buttons according to a user's taste.

| MODE  | SET TIME | TIME    | RF LEVEL | ICE LEVEL | MODE    | SET TIME | TIME    | RF | ICE |
|-------|----------|---------|----------|-----------|---------|----------|---------|----|-----|
| BASIC | 25 min   | 10 min. | 8        | 0         | MANUAL  | 20 min.  | 20 min. | 8  | 2   |
|       |          | 5 min.  | 8        | 3         |         |          | 5 min.  | 7  | 0   |
|       |          | 10 min. | 8        | 0         |         |          | 5 min.  | 7  | 3   |
| LONG1 | 40 min.  | 10 min. | 8        | 0         | SMALL1  | 25 min.  | 5 min.  | 7  | 0   |
|       |          | 5 min.  | 8        | 3         |         |          | 5 min.  | 7  | 3   |
|       |          | 10 min. | 8        | 0         |         |          | 5 min.  | 6  | 0   |
|       |          | 5 min.  | 8        | 3         | SMALL1  | 9 min.   | 3 min.  | 6  | 2   |
|       |          | 10 min. | 7        | 0         |         |          | 3 min.  | 6  | 0   |
|       |          | 5 min.  | 7        | 3         |         |          | 3 min.  | 6  | 2   |
| LONG2 | 40 min   | 15 min. | 7        | 0         | EASY(F) | 10 min.  | 10 min. | 8  | 3   |
|       |          | 10 min. | 7        | 3         |         |          |         |    |     |
|       |          | 15 min. | 7        | 0         |         |          |         |    |     |

